


Digital Microscope VM Slide Explorer 4.0



Fully integrated into Pathology Information Systems and LIMS

- Integration into your LIMS system: Nexus/PATHOLOGIE, PCS/Athena, ...
- Digital slides to a case open automatically, without any click
- Imaging results are transferred to the system automatically



Automatic Quantification of immune stained slides

- Automatic quantification of Ki67, ER PR, CD3, CD15 and more
- Evaluation is performed simultaneously in the background while you work
- Results as image, Excel table or export into your pathology/LIMS system



VMscope GmbH
Campus Charité Mitte
Charitéplatz 1
10117 Berlin, Germany

Phone:
+49 (0)30 450 536188
Fax:
+49 (0)30 450 536910

info@vmscope.com
www.vmscope.com

Brochure: VM Slide Explorer 4
As of: 07th April 2016


© VMscope GmbH Berlin

Digital Microscope VM Slide Explorer 4.0

Tablet Control

Navigation with the tablet PC

- Navigation with the mouse, touch screen or attached tablet PC
- Convenient touch control even at the standard desktop
- Compatible with Apple iPad



Maximum Compatibility


All major formats supported

- VMscope .VSF
- 3DHistech Pannoramic .mrxs
- Hamamatsu .ndpi
- Aperio ScanScope .svs
- Olympus WebImage
- Leica .scn
- Zoomify
- JPEG2000 .jp2
- BigTiff .tif

Multi-purpose Application

For the clinic network or locally

- Evaluation of cases including quantification of tumour markers
- Presentations at medical or clinical conferences using a laptop
- Case collection for routine or studies
- Remote access to cases via networks



VMscope GmbH
Campus Charité Mitte
Charitéplatz 1
10117 Berlin, Germany

Phone:
+49 (0)30 450 536188
Fax:
+49 (0)30 450 536910

info@vmscope.com
www.vmscope.com